

Protecting the Present, *Building for the Future*



Why InSiteOne is Healthcare's Storage Service Provider



- **Full range of data management services:**
 - **Standards based architecture for all stored objects**
 - **One-time storage fee per examination**
 - **A service provider focused only on storage technology**

- **The largest secure online archive for all DICOM objects – radiology, cardiology, mammography**

Internet DICOM Express

The most technologically advanced digital image archive, distribution and management system available.

The Most Powerful Archive

InSiteOne Product Lineup

InDex OnLine

InDex NearLine

InDex Recovery

InDex Web Option

InDex Link


InDex Recovery Plus

InDex® Archive Service



- **An Intelligent DICOM archive service managing onsite spinning archive storage**
- **Up to 7 year onsite capacity to fit workflow requirements !!**
- **All studies stored for indefinitely with verification of all data in archive every sixty days**
- **Includes the disaster recovery solution**
- **Data Migration services without capital expense available**

InDex Online with RIS, Modalities & PACS



Copyright © 2005 InSite One, Inc.
All Rights Reserved

InSiteOne Value Proposition



■ Performance

- Our only business is enterprise storage and data delivery

■ Accessibility

- We support delivery of data 24/7/365

■ Compatibility

- Open architecture standards to accept storage for future clinical systems

InSiteOne Value Proposition



- **Security**
 - We monitor and maintain the highest levels for your data storage architecture

- **Scalability**
 - We take limits out of the question

- **Affordability**
 - We offer an immediate managed storage solution with no capital investment

InSiteOne: The Industry Leader



Protecting the present.....

- **More than 850 million images stored in archive**
- **More than 11 million studies stored**
- **Endorsed by Hospital Associations and GPO's**
- **Partnerships with leading PACS providers**
- **Over 350 customers, over 400 sites**

InSiteOne: The Industry Leader



Building for the future.....

- **InDex 4.0 product upgrade**
- **Dramatic increase in functionality of our products**
- **Increasing PACS provider Partnerships**

The Need for InDex



- **Fast access and remote access to images are modern requirements**
- **Obsolescence and lack of scalability are causing healthcare facilities to search for alternatives to traditional onsite archiving**
- **Regulatory Compliance is driving up the cost and complexity of storing digital medical images**
 - **Storage requirements may be decades on some studies (pediatrics, mammography, OSHA)**
 - **HIPAA requirements include redundancy and disaster recovery**

InDex™ is... Obsolescence Free



- **DICOM-3 standard**
- **Vendor independent**
- **Software and Hardware is continuously improved and upgraded**
- **RAID onsite is InSiteOne provided, maintained, and expanded**

InDex™ is... Scalable



- **Architecture provides for multiple servers per client for virtually unlimited scalability**
- **Multiple servers can act as single Virtual Archive Manager to PACS**
- **Virtually unlimited capacity in data centers available on demand**

InDex™ is... Permanent



- **Archive media is currently DVD**
- **Every disc is error free**
- **Permanent media - cannot be written over or modified**
- **Impervious to hacker attacks, software glitches or viruses**
- **Optical media is more reliable; rated life 100 years**

InDex™ is... High Performance



- **Performance of the archive is dependent on two main factors**
 - **Software and hardware running the archive**
(Archive Manager)
 - **Media which houses the images**
(Spinning disk)

InDex™ is... High Performance



- **InDex uses onsite Linux-based servers capable of processing 1,000,000 studies annually per server**
- **InDex uses RAID on each proxy server for immediate response to image requests**
- **Two levels of intelligent caching minimizes remote retrieval**
- **Lossless compression (up to 3:1) increases throughput**

InDex™ is...




**InDex is the answer for
medical storage, short-
term and long term.**

Thank You!!

Protecting the Present, *Building for the Future*



InDex Online with RIS, Modalities & PACS



Copyright © 2005 InSite One, Inc.
All Rights Reserved