

(800) 554-5090 or (303) 423-7572 or Fax (303) 940-7165

Summit Doppler Systems Inc., 4680 Table Mountain Dr #150, Golden CO 80403 USA - www.summitdoppler.systems

Main Units - The price of each Doppler system includes one probe of your choice.

Probe Choices: See appropriate probe table for probe choice - please specify probe choice when ordering.

Model	Description	Description					Price
		Interchangeable Probes	Bi-Directional Printer	Rechargeable Batteries	Display	Audio Recorder	
L150	Non-Display Hand-Held Doppler	✓					\$545.00
L150R	Non-Display Hand-Held Doppler with Recharger	✓		✓			\$595.00
L150A	Non-Display Hand-Held Doppler with Audio Recorder	✓				✓	\$645.00
L250	Display Hand-Held Doppler	✓			✓		\$745.00
L250R	Display Hand-Held Doppler with Recharger	✓		✓	✓		\$795.00
L250AR	Display Hand-Held Doppler with Recharger and Audio Recorder	✓		✓	✓	✓	\$895.00
L350R	Table-Top Display Doppler with Recharger	✓		✓	✓		\$995.00
L250AB	Reimbursable Vascular Doppler with Portable Printer and Bi-Directional Probe, 2 Cuffs and Aneroid	✓	✓		✓		\$2,190.00
L250AC	Reimbursable Vascular Doppler with Portable Printer and Bi-Directional Probe, 4 Cuffs and Aneroid	✓	✓		✓		\$2,280.00
L450VA	Vista ABI - PVR, 8MHz Probe, Auto Cuff Inflation/Deflation, Stand & 4 cuffs	✓	✓	✓	✓		\$4,995.00
L500VA	Vista AVS - PVR, 8 MHz Probe, PPG, Auto Cuff Inflation/Deflation, Individual Site Mode, PC Software, Stand & 9 Cuffs	✓	✓	✓	✓		\$6,990.00
C250CB	Obstetrical Combination (L250 Display with 2 & 3 MHz Probes)	✓			✓		\$920.00
C250CV	Vascular Combination (L250 Display with 5 & 8 MHz Probes)	✓			✓		\$920.00
C250CC	Obstetrical/Vascular Combo. (L250 Display with 3 & 8 MHz Probes)	✓			✓		\$920.00
C250CP	Obstetrical Combo. (L250 Display with EchoHeart, 2 & 3 MHz Probes)	✓			✓		\$1,200.00
C250CM	Obstetrical Combo. (L250 Display with EchoHeart & 3 MHz Probes)	✓			✓		\$1,050.00

**Non-Directional Probe Table** for L150,L150R,L150A,L250,L250R,L250AR,L350R

Additional Probe Prices

Model	Description		List Price
SD2	2MHz-Obstetrical Probe	Late Term Pregnancy/Obese Patients	\$195.00
SD3	3MHz-Obstetrical Probe (Best general purpose obstetrical use)	Early Term Pregnancy/General Use	\$195.00
SDW	2MHz-Obstetrical Waterproof Probe	Water Applications	\$235.00
	Note: The waterproof can be substituted for OB probes (SD2 & SD3) for an additional charge of:		\$40.00
SD4	4MHz-Vascular Probe	Ultra Wide Beam	\$195.00
SD5	5MHz-Vascular Probe	Pencil Style-Deep Vessels	\$195.00
SD8	8MHz-Vascular Probe (Best general purpose vascular use)	Pencil Style-Surface Vessels	\$195.00
SDE	EchoHeart 5MHz Transvaginal Obstetrical Probe	Fetal Heart detection as early as 6 wks.	\$395.00
	Note: The EchoHeart can be substituted for OB probes (SD2,SD3 & SDW) for an additional charge of:		\$200.00

**Bi-Directional Probe Table** for L250AB,L250AC, L450VA & L500VA

Additional Probe Prices

SD5B	5MHz-Bi-Directional Vascular Probe for L250AB & L250AC	Pencil Style-Deep Vessels	\$295.00
SD8B	8MHz-Bi-Directional Vascular Probe for L250AB & L250AC	Pencil Style-Surface Vessels	\$295.00
SD5V	5MHz-Bi-Directional Vascular Probe for L450VA & L500VA	Pencil Style-Deep Vessels	\$395.00
SD8V	8MHz-Bi-Directional Vascular Probe for L450VA & L500VA	Pencil Style-Surface Vessels	\$395.00

**Additional Accessories**

Part No.	Description		List Price
L450PG	PPG Package - includes PPG Probe & 1.9 cm Cuff	For L450VA - Vista ABI	\$495.00
L450SF	Software Package - includes PC Software & Interface	For L450VA - Vista ABI	\$795.00
K150	ABI Kit-with Aneroid, 2 Cuffs, Forms, Chart, ABI Video	Not For L250AB or L250AC	\$195.00
K160	ABI 2 Cuff Package - 10 or 12cm & Aneroid	Vascular Uses - ABI Exams	\$175.00
K170	ABI 4 Cuff Package - 10 or 12cm & Aneroid	Vascular Uses - ABI Exams	\$245.00
K180	ABI Printer Label Paper - 5 Rolls (30 Waveforms Per Roll)	For L250AB, L250AC, L450VA, L500VA	\$34.00
K200	Stand with Storage Basket - For Hand-Held Systems	For L150's/L250's	\$245.00
K220	Stand with Storage Basket - For Table-Top Systems	For L350R	\$245.00
K230	Stand without Storage Basket - For Hand-Held Systems	For L150's/L250's	\$195.00
K240	Stand without Storage Basket - For Table-Top Systems	For L350R	\$195.00
K260	Carrying Case for all Hand-Held Dopplers	For L150's/L250's	\$19.95
CBL0008	Security Locking Cable	For L150's/L250's/L350's	\$49.95
G150	Ultrasound Gel 60 gm	All Dopplers	\$2.95

Service 1-3 day Guarantee Repair Turn-Around Time

Part No.	Description	MAXIMUM
SVMAIN	Service-Main Unit After One Year Warranty Period	\$149.00
SVPROB	Service-Probe After One Year Warranty Period	\$79.00