

# Mammography Cape Warmers



- Provide your mammography patients with the relaxing comfort of a warm cape
- Made of 600 Denier Nylon with Dupont Thermoloft insulation
- Thermostat controlled heating elements (top and bottom)
- Warmer interior maintains a safe, consistent temperature
- Heating elements operating temperature, 130 degrees
- 24/7 operation recommended
- Operates on 110 volts AC
- 4.5 Amps power consumption
- Weight: 3 - 4.3 lbs
- Available in three sizes

MRXW-1LS-AC	Mammography Cape Warmer 18x19x5.5 with AC Power Supply
MRXW-XWLS-AC	Mammography CapeWarmer 19x21x8.5 with AC Power Supply
MRXW-2LS-AC	Mammography Cape Warmer 18x19x12 with AC Power Supply

