

LAB COATS Premium Disposable

Medgluv disposable premium lab coats are engineered to optimize comfort and barrier protection. Lab coats are rigorously tested and manufactured in an ISO Certified facility under stringent process controls to ensure that each product meets exacting quality standards and performs to specification

Characteristics and Overview:


- Knee-length provides ample body coverage.
- Knit collar, knit cuffs and snap front
- Elastic wrists provide a secure fit and reduce the risk of tears and snags
- Made with protective multilayer nonwoven material that is fluid resistant, breathable and comfortable
- Triple layer Non-Conductive Material resists moisture dirt and spill.



Mask Specifications:

Brand: Disposable Premium Lab Coats
 Material: Spunbond non-woven material
 Color: Blue and White
 Item Number: MGLC10 / MGLC20
 Packaging: 10pcs/bag, 3bags/case

WITHOUT POCKETS			WITH 3 POCKETS		
Item #	Size	Color	Item #	Size	Color
MGLC101	Small	White	MGLC201	Small	White
MGLC102	Medium	White	MGLC202	Medium	White
MGLC103	Large	White	MGLC203	Large	White
MGLC104	X-Large	White	MGLC204	X-Large	White
MGLC105	2X-Large	White	MGLC205	2X-Large	White
MGLC106	Small	Blue	MGLC206	Small	Blue
MGLC107	Medium	Blue	MGLC207	Medium	Blue
MGLC108	Large	Blue	MGLC208	Large	Blue
MGLC109	X-Large	Blue	MGLC209	X-Large	Blue
MGLC110	2X-Large	Blue	MGLC210	2X-Large	Blue



Single use only
Do not reuse.



Latex-Free



Non-Sterile

