







VS30 AND VS400 COMPARISON SHEET

ITEM		VS30	VS400
Appearance	Main unit		
	Charging cradle	No	
	Charging adapter		
	Stand		
		Both fixed stand and mobile stand of VS30 are compatible with VS400.	
Function	Color	3 colors, Cyan, White, Red	5 colors, Cyan, White, Red, Blue, Green.
	Brightness	3 degrees	3 degrees
	Size	2 size	3 size
	Inverse	Inverse mode is available	Inverse mode is available
	Enhance mode	No enhance mode	Enhance mode is available

ITEM		VS30	VS400
Specification	Vein Detecting Light Source	Near-infrared light	
	Infrared Light Wavelength	Approximately 850 nm	
	Visible Light Wavelength	Approximately 530 nm	
	Light Radiation	No	
	Battery	Lithium-ion Battery, 3.6V, 4300mAh	
	Continuous Run Time	4 hours	5 hours
	Battery Charge Time (Full)	3 hours	
	Power Supply Input	100-240VAV, 50-60Hz, 0.4A	
	Battery Charge Input	5VDC, 2A	
	Weight	460g	480g
	Dimension (L*W*H)	240mmx60mmx60mm	
	Focusing Distance	Approximately 20cm	