

OysterSkin 250 Nitrile Exam Gloves

Powder-Free Non-Sterile Greay Nitrile Exam Gloves

OysterSkin Nitrile Exam Gloves provide tactile sensitivity and protects hands from potentially infectious materials. New nitrile compounding technology allows for a thinner glove with incredible elasticity. They offer a soft and comfortable feel with outstanding tactile sensitivity.

Packaged in 250 gloves per box.

Devotes less space to glove storage, less time handling cartons.

Characteristics and Overview:

- Finger-tip textured enable grip even when wet.
- Virtually no pressure on the hands and reduces fatigue after extended wear.
- Offer excellent sensitivity for touch perception and dependable protection.
- 60% saving in the storage space, fewer boxes to discard and less waste.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.
- Meets and exceeds NIOSH recommendations for Fentanyl.



Item No.	Description	Size	Packaging
MG5301	OysterSkin Nitrile PF Grey Glove	X-Small	250 eaches/box 10 boxes/case
MG5302	OysterSkin Nitrile PF Grey Glove	Small	250 eaches/box 10 boxes/case
MG5303	OysterSkin Nitrile PF Grey Glove	Medium	250 eaches/box 10 boxes/case
MG5304	OysterSkin Nitrile PF Grey Glove	Large	250 eaches/box 10 boxes/case
MG5305	OysterSkin Nitrile PF Grey Glove	X-Large	230 eaches/box 10 boxes/case

Glove Specifications:

Brand	OysterSkin
Material	Nitrile
Color	Grey
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	XS, SM, MD, LG and XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Textured Finger-tips
Length (inches)	9.5 in
Palm Thickness (Avg)	3.1 mil
Finger Thickness (Avg)	3.5 mil
Tensile Strength (Mpa)	35 / 40 Mpa
Ultimate Elongation (%)	560 / 550
Pin Hole (AQL)	1.0

Summary of Chemotherapy Testing Reports – ASTM D6978-05	
Bleomycin Sulfate 15mg/mL, >240 mins.	Idarubicin 1mg/mL, >240 mins.
Busulfan 6mg/mL, >240 mins.	Ifosfamide 50mg/mL, >240 mins.
Carboplatin (Paraplatin) 10mg/mL, >240 mins.	Irinotecan 20mg/mL, >240 mins.
Carmustine (BCNU) 3.3mg/mL, 17.1 mins.	Mechlorethamine 1mg/mL, >240 mins.
Cisplatin 1.0mg/mL, >240 mins.	Melphalan 5mg/mL, >240 mins.
Cyclophosphamide (Cytosan) 20mg/mL, >240 mins.	Methotrexate 25mg/mL, >240 mins.
Cytarabine 100 mg/mL, >240 mins.	Mitomycin C 0.5mg/mL, >240 mins.
Dacarbazine (DTIC) 10mg/mL, >240 mins.	Mitoxantrone 2mg/mL, >240 mins.
Daunorubicin 5mg/mL, >240 mins.	Paclitaxel (Taxol) 6mg/mL, >240 mins.
Docetaxel 10mg/mL, >240 mins.	Rituximab 10mg/mL, >240 mins.
Doxorubicin Hydrochloride 2mg/mL, >240 mins.	Thiotepa 10mg/mL, 27.9 mins.
Epirubicin (Ellence) 2mg/mL, >240 mins.	Trisenox 1mg/mL, >240 mins.
Etoposide (Toposar) 20mg/mL, >240 mins.	Vincristine Sulfate 1mg/mL, >240 mins.
Fludarabine 25mg/mL, >240 mins.	
Fluorouracil 50mg/mL, >240 mins.	Fentanyl tested as follows:
Gemcitabine 38mg/mL, >240 mins.	Fentanyl citrate 100mcg/2mL, >240 mins.
WARNING: Not for use with Carmustine or Thiotepa.	

Quality Standards:

Manufactured in accordance with Quality System ISO 9001.

Exceeds and meets current ASTM Standards: - D6319, D5151, D1671 & D6978.

This product is made from Synthetic Nitrile Polymer and contains no natural latex rubber.

* When used in accordance to the guidelines set forth in the USP 800 chapter.

Tested chemotherapy drugs available upon request.