

# The PhysioPRO II



**The PhysioPRO II a portable, clinical grade shockwave therapy system that induces an intracellular and extracellular reactions to assist in the healing and recovery of soft tissue.**

**Shockwave technology theory works by a deep penetrating focused wave that activate cells, driving the formation of new cells that replace damaged tissues.**

**Shockwave therapy accelerates blood flow to the treated area and may induce vasodilation and angiogenesis restoring and improving musculoskeletal tissue.**

# The PhysioPRO II

## The PhysioPRO II: Portable Wave Therapy for Pain Management

### Pre-Programmed Therapy Protocols for:

- Radial And Ulnar Epicondylitis
- Tendinitis Of The Shoulder
- Acute Calcific Tendinitis Of The Shoulder
- Post Muscular Injury
- Patellar Tendinopathy
- Achillodynia
- Plantar Fasciitis
- Heel Spurs
- Myofascial Trigger Point Therapy - Neck
- Myofascial Trigger Point Therapy - Back
- Trochanteric Bursitis
- Periostitis / Shin Splints
- Dupuytren's Disease
- Thumb Basal Joint Arthritis / Rhizarthritis

### System Features & Benefits:

- Controlled manually or remotely by app for phone or tablet.
- "Plug & play" system
- Compact and Portable
- Foot Pedal Not Required
- Multi-preset treatment protocols
- User friendly Interface
- Low Cost of Ownership

### The PhysioPRO II TECHNICAL SPECIFICATIONS:

<b>System Type</b>	Main Unit with Therapy handle	<b>Protocols</b>	More than 25 illustrated preset treatment recommendations
<b>Technology</b>	Compressor free ballistic radial pulse therapy-system with electromagnetic generator	<b>Controls</b>	Phone or tablet app for all software operations
<b>Power Levels / Energy</b>	60 to 185 mJ (equivalent to 1-5 bar)	<b>Applicators</b>	6/15/25mm
<b>Frequency</b>	1-22 Hz	<b>Life time</b>	Minimum warranty of 4,000,000 shocks per handpiece Minimum warranty of 150,000 shocks per applicator head
<b>Modes</b>	4 Burst modes (continuous / 4 / 8 / 12 pulses)	<b>Dimensions</b>	12 x 9.5 x 5 Inches (L/W/H)
<b>Programs</b>	7 preset programs, adjustable	<b>Weight</b>	2.07kg (controller unit only)
<b>Package</b>	Aluminum Carrying Case	<b>Warranty</b>	3 years for controller unit and 1 year for accessories

### The PhysioPRO II STANDARD ACCESSORIES:

ORDERING CODE	DESCRIPTION	QUANTITY / PC.	ORDERING CODE	DESCRIPTION	QUANTITY / PC.
OCE001010101Z	Therapy Handle	1	OCE00102040001	Power Line	1
OCE001020102	Foot Switch	1	OCE001010101Z-04	Handle Pedestal	1
OCE001010101Z-07	6mm Therapeutic Head	1	OCE00101010405	Stylus Pen	1
OCE001010101Z-08	15mm Therapeutic Head	1	OCE001010101Z-06	Silicone Cap	10
OCE001010101Z-09	25mm Therapeutic Head	1	OCE00103040101	Aluminum Alloy Suitcase	1

Florida Medical Sales

4320 A1A South, Suite 2, Saint Augustine, FL 32080

admin@floridamedicalsales.com | Tel: 888.269.1810