



Checkme O2 Pro

Smart wrist oximeter with memory function

Compatible for Overnight Studies with PC Download and Report



The Viatom Checkme O2 Pro is the next generation pulse oximeter for overnight monitoring of SpO2 and pulse rate. Compatible with PC-Download, for complete report printout. The pulse oximeter uses a rechargeable battery and uses a revolutionary patented 'ring' pulse oximeter, for the best patient experience we have encountered.


Part Number : VIAPRO-01

Includes the following:

Wrist oximetry monitor with ring sensor and vibration
Sleep mode, Monitor mode, Bluetooth, APP, PC software

- Four night recording memory
- Easy-to-clean, reusable ring sensor or other sensors available
- Rechargeable battery
- Compatible with PC-Download
- Soft, comfortable silicone wrist strap
- Nail polish? No problem! Sensor fits around base of thumb away from nails

The Viatom Checkme O2 Pro is the next generation wrist-worn pulse oximeter. Designed from the bottom up for the best patient experience, with a simple wrist strap, single button, and download cable, this device is easy to use and to clean. The reusable ring sensor is the device's secret weapon- incredibly comfortable, your patient will probably forget they are wearing a monitoring device at all.



Small, light, with 24 hours battery and room for four recordings. Bluetooth Compatible with smartphones for real time display and recordings.

Patented Ring sensor (included) never comes off the patient's thumb, is comfortable and soft, and should lead to a better night's sleep than previous devices.

Technical specification

Recorded channels	SpO2 , pulse rate, motion
Dimensions	44 x 25 x 15 mm
Weight including batteries.....	12 g (main unit)
Power	Internal rechargeable batteries. 24 hours expected for typical recording parameters
Sensors.....	Ring sensor included as standard.
Data storage	Internal memory holds four nights' recordings
Setup & signal check	Via OLED screen
Software	Download from website. Custom reports and SPO2 threshold configurable
CE.....	Viatom holds CE 0197
Warranty.....	12 months on main unit, 6 months on probes