

Elite XR Series



Elite XR Comfort Series

- VisionTec™ lenses feature RTC™ technology (Radiation Transmission Control™) Attenuation: 95% radiation reduction at 125 KVP direct beam.
- Anti-fog coating prevents fogging and makes cleaning easier
- Multi layered anti-reflective coating increase visual acuity and performance
- HLTO™ (High Light Transmission Optics™) and the VisionTec™ lens offer superior optical clarity, one of the highest visible light transmittance in the industry allowing you to see better in the procedure environment while maintaining uncompromised radiation protection
- Rx'able (Prescription Available - single vision and bifocal)
- Lead equivalent: 0.75 (+/- .05) Pb glass

Frame:

- Ergonomic frame design
- Durable TR90 frame material
- Available in black, pink, and dark grey

Accessory Bonus:

- Micro fiber cleaning cloth
- Protective case
- Anti-Reflection Spray Cleaner

All Vision Scientific™ products meet ANSI Z87.1 Safety Standards and is backed by a one year limited warranty against manufacturer's defect from date of purchase.

Pricing does not include shipping.



Elite XR Fitover Series



Elite XR Fitover Series

- VisionTec™ lenses feature RTC™ technology (Radiation Transmission Control™) Attenuation: 95% radiation reduction at 125 KVP direct beam.
- Anti-fog coating prevents fogging and makes cleaning easier
- Multi layered anti-reflective coating increase visual acuity and performance
- HLTO™ (High Light Transmission Optics™) and the VisionTec™ lens offer superior optical clarity, one of the highest visible light transmittance in the industry allowing you to see better in the procedure environment while maintaining uncompromised radiation protection
- Rx'able (Prescription Available - single vision and bifocal)
- Lead equivalent: 0.75 (+/- .05) Pb glass

Frame:

- Fitover frame design
- Durable TR90 frame material
- Extra wide temple allow for a comfortable fit with your existing eyewear
- Available in black, demi, and dark grey

Accessory Bonus:

- Micro fiber cleaning cloth
- Protective case
- Anti-Reflective Spray Cleaner

All Vision Scientific™ products meet ANSI Z87.1 Safety Standards and is backed by a one year limited warranty against manufacturer's defect from date of purchase.

Pricing does not include shipping.



Elite XR Metal Series



Elite XR Metal Series

- VisionTec™ lenses feature RTC™ technology (Radiation Transmission Control™) Attenuation: 95% radiation reduction at 125 KVP direct beam.
- Anti-fog coating prevents fogging and makes cleaning easier
- Multi layered anti-reflective coating increase visual acuity and performance
- HLTO™ (High Light Transmission Optics™) and the VisionTec™ lens offer superior optical clarity, one of the highest visible light transmittance in the industry allowing you to see better in the procedure environment while maintaining uncompromised radiation protection
- Rx'able (Prescription Available - single vision and bifocal)
- Lead equivalent: 0.75 (+/- .05) Pb glass

Frame:

- Ergonomic frame design
- Perfect Fit System™ with flex fix spring hinges increase comfort
- Nickel Free Coating - hypoallergenic
- Adjustable silicone nose pad system
- High grade ultra-light Aluminum Alloy (Atomic #13) on XR Series frames.
- Monel Alloy - durability and strength.
- Available in gunmetal and matte brown

Accessory Bonus:

- Micro fiber cleaning cloth
- Protective case
- Anti-Reflective Spray Cleaner

All Vision Scientific™ products meet ANSI Z87.1 Safety Standards and is backed by a one year limited warranty against manufacturer's defect from date of purchase.

Pricing does not include shipping.