

# KwikVac

## SMOKE FREE EVACUATION

### Procedure Operating Modes

The desired procedure can be set simply by touching the corresponding image. Automatically customizing suction default level based on desired surgical procedure.

### Longer lasting filter life (up to 50 Hours)

Can save on running costs.



### Touch screen

The touch screen panel can be operated while wearing gloves.

### Multi-Port design

Three different sized connection ports are designed to create an air-tight bond with the desired port. The port cover cap eliminates any possible contamination with the glove when connecting the hose to the unit.

## Technical Specification

Item	Specification
Noise Levels (Nominal Values)	maximum power < 55 dBA
Dimensions (H x W x D)	6 in x 11 in x 15.5 in (15.2 cm x 27.9 cm x 39.4 cm)
Flow Rates	35 cfm (991 lpm) using 7/8 in (22 mm) tubing
Weight	With filter installed: 12.0 lbs (5.0 kg) Without filter installed: 9.5 lbs (4.0 kg)
Warranty	Three (3) Years Full Warranty
Certifications	CE, FDA
Multi Port	1/4 in (6.4mm), 3/8 in (9.5mm), 7/8 in (22mm)

Model	Product Description	Included	Packaging Info.
KWV-M	KwikVac Electro Surgical Evacuator (Without filter, Includes one footswitch)	-	1 / each
KWV-F	KwikVac Filter, Multi ports design with ULPA filter(ULPA particle Filter with 99,999+% efficiency)	Optional	4 / Case
KWV-DT	KwikVac Detecto, Automatic Activation Device	Optional	1 / Box
KWV-LT	KwikVac Liquid Trap, 7/8 in (22 mm) and 1/4 in (6.4 mm) tubing connection ports	Optional	20 / Case

## Accessories



Foot-switch (Included)



Detecto Activation Device (Optional)




Fluid Trap (Optional)

# KwikVac Smoke Evacuation Pencils

Extendable Smoke Evacuation Pencil

One Touch Button

Extendable Type




Model	Product Description	Type	Cable Length	Packaging Info.
EXTL	Extendable, 2.5" Non-Stick Ceramic Coated Blade, Holster	Swivel	10' Cable	1 pc/ Pack 20pcs/ Box
EXTL-5	Extendable, 2.5" Non-Stick Ceramic Coated Blade, Holster	Swivel	16' Cable	
EXTLI	Extendable, 2.5" Insulated Non-Stick Ceramic Coated Blade, Holster	Swivel	10' Cable	
EXTLI-5	Extendable, 2.5" Insulated Non-Stick Ceramic Coated Blade, Holster	Swivel	16' Cable	

## Regular Smoke Evacuation Pencil

**COST EFFECTIVE and SIMPLE TO USE**



**Extension of the suction tip** allows maximum smoke evacuation



Model	Product Description	Type	Cable Length	Packaging Info.
REGL	Regular, 2.5" Non-Stick Ceramic Coated Blade, Holster	Swivel	10' Cable	1 pc/ Pack 20pcs/ Box
REGL-5	Regular, 2.5" Non-Stick Ceramic Coated Blade, Holster	Swivel	16' Cable	
REGLI	Regular, 2.5" Insulated Non-Stick Ceramic Coated Blade, Holster	Swivel	10' Cable	
REGLI-5	Regular, 2.5" Insulated Non-Stick Ceramic Coated Blade, Holster	Swivel	16' Cable	

### Non-Stick Blade – Easy to Clean


Ceramic coated blade is specially designed for easy cleaning after cutting and coagulation. The blade is not flammable, inhibits noxious gasses, bio-compatible to the skin and extremely durable when using higher wattages.



Before use



After use in 50W




Easy cleaning with one swipe

# WelFit

## Surgical Smoke Evacuation Adapter

Extension of the suction tip allows maximum smoke evacuation



Model	Product Description	Tubing Length	Packaging Info.
WELF	Ideal for surgical smoke evacuator applications, Extension	3/8in X 10' Tubing	1 pc / Pack
WELF-5	Ideal for surgical smoke evacuator applications, Extension	3/8in X 16' Tubing	20 pcs/ Box

# **MediSourceNet**

1303 W Valencia Dr #243

Fullerton CA 92833

Email : [Info@Medisourcenet.com](mailto:Info@Medisourcenet.com)

Tel : 657-243-2848