

# ELECTRODE PRODUCT GUIDE

See reverse side for proper cable/connector requirements.

Product Name & Number	Differentiation Features	Location Used	Defibrillator
<b>stat•padz</b> (8900-4003) 	<ul style="list-style-type: none"> <li>Green pressure sponge - Promotes better skin coupling.</li> <li>Pre-Connect - Saves critical application time.</li> <li>Speed pack - Organizes wires to eliminate tangles and fumbling.</li> <li>Different pad shapes - Visual cues for proper pad placement.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series®</li> <li>E Series®</li> <li>R Series®</li> </ul>
<b>stat•padz II</b> (8900-0802-01) 	<ul style="list-style-type: none"> <li>Connects directly to AEDs and ALS devices.</li> <li>See <b>stat•padz®</b> for other benefits.</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>AED Plus</li> <li>AED Pro</li> </ul>
<b>CPR stat•padz</b> (8900-0400) 	<ul style="list-style-type: none"> <li>Designed for ALS market, two-piece design ensures all therapies can be delivered.</li> <li>Provides defibrillator with Real CPR Help®.</li> <li>Accommodates AP and AA pad placement.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>AED Plus</li> <li>AED Pro</li> </ul>
<b>CPR-D•padz</b> (8900-0800-01) 	<ul style="list-style-type: none"> <li>Provides defibrillator with Real CPR Help.</li> <li>One-piece design simplifies proper positioning on the chest.</li> <li>One size fits all, with ability to expand for larger chest.</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> <li>Public access (i.e. lobby)</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>AED Plus</li> <li>AED Pro</li> </ul>
<b>pro•padz Radiolucent Solid Gel</b> (8900-4005) 	<ul style="list-style-type: none"> <li>For use in Cath and EP labs for procedures - Virtually invisible to X-ray.</li> <li>Pre-Connect - Saves critical application time.</li> <li>Speed pack - Organizes wires to eliminate tangles and fumbling.</li> <li>Different pad shapes - Visual cues for proper pad placement.</li> </ul>	<ul style="list-style-type: none"> <li>Cardiac catheterization lab</li> <li>Electrophysiology labs</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pro•padz Radiolucent Liquid Gel</b> (8900-2105-01) 	<ul style="list-style-type: none"> <li>For use in Cath and EP labs for procedure.</li> <li>Liquid gel provides optimal skin coupling for no dermal damage.</li> <li>Virtually invisible to X-ray.</li> </ul>	<ul style="list-style-type: none"> <li>Cardiac catheterization lab</li> <li>Electrophysiology labs</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pro•padz Sterile</b> (8900-4012) 	<ul style="list-style-type: none"> <li>For use in operating rooms.</li> <li>Helps maintain sterile barrier.</li> <li>Packaged in two pouches.</li> </ul>	<ul style="list-style-type: none"> <li>Operating rooms</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pro•padz Cardiology Specialty</b> (8900-2100-01) 	<ul style="list-style-type: none"> <li>No-slump tray.</li> <li>Liquid gel provides optimal skin coupling for no dermal damage.</li> </ul>	<ul style="list-style-type: none"> <li>Cardiac catheterization lab</li> <li>Electrophysiology labs</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pro•padz Biphasic</b> (8900-2303-01) 	<ul style="list-style-type: none"> <li>Used in Cath and EP labs for procedures.</li> <li>Highest conversion rate for cardioversions with low energy biphasic waveform.</li> <li>Liquid gel provides optimal skin coupling for no dermal damage.</li> </ul>	<ul style="list-style-type: none"> <li>Cardiac catheterization labs</li> <li>Electrophysiology labs</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pedi•padz</b> (8900-2065) 	<ul style="list-style-type: none"> <li>Meets all pediatric ALS criteria.</li> <li>Liquid gel provides optimal skin coupling for no dermal damage.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pedi•padz Multi-Function Solid Gel</b> (8900-3000-01) 	<ul style="list-style-type: none"> <li>Meets all pediatric ALS criteria.</li> <li>Solid gel provides ease for storing.</li> <li>Pre-connect saves critical application time.</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pedi•padz Radiolucent Solid Gel</b> (8900-1005-01) 	<ul style="list-style-type: none"> <li>Meets all pediatric ALS criteria.</li> <li>For use during procedures.</li> <li>Virtually invisible to X-ray.</li> </ul>	<ul style="list-style-type: none"> <li>Cardiac catheterization lab</li> <li>Electrophysiology labs</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul>
<b>pedi•padz II</b> (8900-0810-01) 	<ul style="list-style-type: none"> <li>Pediatric intelligence function only available on ZOLL® AEDs with <b>pedi•padz® II</b>.</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> <li>Public access (i.e. lobby)</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>AED Plus</li> <li>AED Pro</li> </ul>
<b>OneStep Basic</b> (8900-0211-01) 	<ul style="list-style-type: none"> <li>Ready Pack packaging.</li> <li>Code-Ready® connector.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> </ul>
<b>OneStep Pacing</b> (8900-0212-01) 	<ul style="list-style-type: none"> <li>OneStep® Basic features plus:</li> <li>Eliminates need for separate ECG leads.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> </ul>
<b>OneStep CPR</b> (8900-0213-01) 	<ul style="list-style-type: none"> <li>OneStep Basic features plus:</li> <li>Enable Real CPR Help with visual and verbal cues.</li> <li>See-Thru CPR®.</li> <li>Anterior - Posterior placement.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> </ul>
<b>One Step CPR (AA)</b> (8900-0217-01) 	<ul style="list-style-type: none"> <li>OneStep Basic features plus:</li> <li>Enable Real CPR Help with visual and verbal cues.</li> <li>See-Thru CPR.</li> <li>Anterior - Anterior placement.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> </ul>
<b>OneStep Complete</b> (8900-0214-01) 	<ul style="list-style-type: none"> <li>All OneStep Basic, Pacing, and CPR capabilities.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> </ul>
<b>OneStep Pediatric</b> (8900-0215-40) 	<ul style="list-style-type: none"> <li>Designed for patients 0 - 8 years of age and less than 55 pounds (25 kg).</li> <li>Automatically adjusts joules setting to pediatric levels.</li> </ul>	<ul style="list-style-type: none"> <li>Med/surge floors</li> <li>ER &amp; ICU</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> </ul>

# ZOLL CPR ELECTRODE GUIDE

See reverse side for individual product details.

**CPR-D • padz**  
(8900-0800-01)

**CPR stat • padz**  
(8900-0400)

**OneStep CPR**  
(8900-0213-01)

**OneStep CPR AA**  
(8900-0217-01)

**OneStep Complete**  
(8900-0214-01)

**OneStep Pediatric**  
(8900-0215-40, 8 per case)  
(8900-0218-40, 1 per case)

**OneStep Pacing**  
(8900-0212-01)



## AEC PLUS



- Connects directly for Real CPR Help® and defibrillation

- Connects directly for Real CPR Help and defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

- Not compatible

## AED PRO



- Connects directly for Real CPR Help and defibrillation

- Connects directly for Real CPR Help and defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

- Not compatible

## M Series



- Requires CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defibrillation

- Requires CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defib/pace

- Requires M Series® CPR Cable for CPR (pt# 8000-0599)
- Connects directly for defibrillation

- Requires M Series CPR Cable for CPR (pt# 8000-0599)
- Connects directly for defibrillation

- Requires M Series CPR Cable for CPR (pt# 8000-0599)
- Connects directly for defibrillation

- Requires M Series CPR Cable for CPR (pt# 8000-0599)
- Connects directly for defibrillation

- Requires M Series CPR Cable for CPR (pt# 8000-0599)
- Connects directly for defib/pace

## E Series



- Requires CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defibrillation

- Requires CPR Connector (pt# 8000-0370) for CPR
- Connects directly for defib/pace

- Not compatible

- Not compatible

- Not compatible

- Not compatible

- Not compatible

## R Series



- Requires R Series® CPR Connector for CPR (pt# 8009-0020)
- Connects directly for defibrillation

- Requires R Series CPR Connector for CPR (pt# 8009-0020)
- Connects directly for defib/pace

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

Function	Defibrillate Only	Multifunction	Multifunction	Multifunction	Multifunction	Multifunction	Multifunction
Construction	1 piece	2 piece	2 piece	2 piece	2 piece	2 piece	2 piece
Pad Placement	AA	AA/AP	AP	AA	AP	AP	AP
Shelf Life	5 yrs	2 yrs	2 yrs	2 yrs	2 yrs	2 yrs	2 yrs
Quantity Sold	Individual	Individual or 8 per case	8 per case	8 per case	8 per case	Individual or 8 per case	8 per case