







Capnoline / Oridion Compatible Sidestream EtCO2 Sampling Line

Condition	Ref. Picture	Connector (Device side)	P/N	Description	Humidity	Patient Type	Compatible With:
Intubated			300.010.002	Intubated EtCO2 sampling line for short-term use	Low	Adult	Microstream® Ext: M3015A, M3015B, MP5: M8105A, M8105AS, VM8: 863066, VM1: 863266Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
Non-Intubated			300.010.001	Nasal Sampling Line for short-term use	Low	Adult	
Non-Intubated			300.010.003	Nasal Sampling Line with O2 tube for short-term use	Low	Adult	