

SCATTERGUARD

DISPOSABLE CAPS



THE PROTECH PROMISE

Protect yourself and your team with Scatterguard caps. ScatterGuard products provide excellent protection for fluoro-lab and operating suite personnel during interventional procedures. They are lightweight, disposable, and easily adjusted to fit all head sizes. ScatterGuard caps are also lead free and environmentally safe for disposal through normal waste procedures.

“The scattered radiation from the patient comprises the main source of radiation dose to staff.” – International Atomic Energy Agency

THE SCATTERGUARD CAP COMES IN 3 PROTECTIVE MODELS:

SG-DCAP-125: 0.125mm LE at 75% Attenuation

SG-DCAP-250: 0.25mm LE at 90% Attenuation

SG-DCAP-375: 0.375mm LE at 95% Attenuation

PROTECH MEDICAL INC