

NitraPro 16

16" Powder-Free Non-Sterile Nitrile Examination Gloves

Excellent for sterilization, decontamination, disinfect, cleaning or handling with chemicals in a hazardous environment. This 16-inch specialty glove is designed to offer exceptional added coverage and protection where there is potential exposure to sharps and chemicals.

Characteristics and Overview:

- Full textured fingertips enhance grip in wet and dry situations.
- Longer 16" extended length
- Tested for puncture resistance and color-coded sizing on box.
- Meets both ASTM F1342 and D6978-05 Standards for use of chemotherapy gloves when handling hazardous drugs*

Item No	Description	Size	Packaging
MG50161	NitraPro Nitrile PF Exam Glove	Small	50 eaches/box 10 boxes/case
MG50162	NitraPro Nitrile PF Exam Glove	Medium	50 eaches/box 10 boxes/case
MG50163	NitraPro Nitrile PF Exam Glove	Large	50 eaches/box 10 boxes/case
MG50164	NitraPro Nitrile PF Exam Glove	X-Large	50 eaches/box 10 boxes/case



Glove Specifications:

Brand	NitraPro
Material	Nitrile
Color	Blue
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	XS, SM, MD, LG, XL and 2XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free
Glove Exterior	Full Textured
Length (inches)	16 in
Palm Thickness (Avg)	7.0 mil
Finger Thickness (Avg)	8.9 mil
Tensile Strength (Mpa)	32 / 30 Mpa
Ultimate Elongation (%)	680 / 620
Pin Hole (AQL)	1.5



Quality Standards:

Manufactured in accordance with Quality System ISO 9001.

Exceeds and meets current ASTM Standards: - D6319, D5151, D1671 & D6978.

This product is made from Synthetic Nitrile Polymer and contains no natural latex rubber.

* When used in accordance to the guidelines set forth in the USP 800 chapter.

Tested chemotherapy drugs available upon request.

