




Introducing SunTech's Equine Blood Pressure and Vital Signs Monitor



The SunTech Vet25E and Vet30E offer all the great features of the Vet25 and Vet30, with the addition of a new blood pressure (BP) algorithm specifically tuned to accurately measure BP in equine patients.





The Vet25E and Vet30E come with cuff sizes #1 through #8

Coming Soon!

Reusable, cleanable cuffs in equine sizes

Specifications	Vet25E	Vet30E
Systolic	40 to 265 mmHg	40 to 265 mmHg
Diastolic	20 to 200 mmHg	20 to 200 mmHg
Mean Arterial Pressure	27 to 222 mmHg	27 to 222 mmHg
Heart Rate	15 to 300 BPM	15 to 300 BPM
Device Weight	2.1 lbs (with battery) 953 grams (with battery)	2.65 lbs (with battery) 1202 grams (with battery)
Dimensions	6.25" x 5" x 5x25" 15.9 cm x 12.7 cm x 13.3 cm	6.25" x 5" x 5x25" 15.9 cm x 12.7 cm x 13.3 cm
Warranty	2 years on monitor 90 days on cuffs/hoses	2 years on monitor 90 days on cuffs/hoses 1 year on AccuVet accessories 6 months on Masimo accessories
Temperature Sensor		YSI 400 Compatible
Temperature Range		0°C to 50°C
SpO ₂ Options		Masimo™ or Accuvet

All SunTech monitors are available in 3 colors:



Flamingo Pink



Peacock Blue



Tree Frog Green