

NitraGrip Pro Nitrile Exam Gloves

12-inch Powder-Free 2-Ply Blue/White pH Examination Gloves

These 2-Ply technology outer Blue and inner White elongated latex-free nitrile gloves offers the unique ability to detects tears and rips by showing the high visibility white inner through the breach.

Characteristics and Overview:

- Designed specially for the Emergency Medical Profession for high risk situation that requires the highest standards in double protection and safety in administering critical and patient care.
- Ideal for EMS, Law Enforcement and Corrections
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM F1342 and D6978-05 Standards for use of chemotherapy gloves when handling hazardous drugs*
- Meets and exceeds the standard emergency glove medical of NFPA 1999 Standards on protective clothing, Emergency Medical Operations, Edition 2018.



| Item No | Description | Size | Packaging |
|---------|----------------------------------|----------|--------------------------------|
| MG50091 | NitraGrip Pro Nitrile Exam Glove | Small | 50 eaches/box 10 boxes/case |
| MG50092 | NitraGrip Pro Nitrile Exam Glove | Medium | 50 eaches/box 10 boxes/case |
| MG50093 | NitraGrip Pro Nitrile Exam Glove | Large | 50 eaches/box 10 boxes/case |
| MG50094 | NitraGrip Pro Nitrile Exam Glove | X-Large | 50 eaches/box 10 boxes/case |
| MG50095 | NitraGrip Pro Nitrile Exam Glove | 2X-Large | 45 eaches/box 10 boxes/case |
| MG50096 | NitraGrip Pro Nitrile Exam Glove | 3X-Large | 45 eaches/box 10 boxes/case |

Glove Specifications:

| | |
|-------------------------|------------------------------|
| Brand | NitraGrip Pro |
| Material | Nitrile |
| Color | Blue (outer) / White (inner) |
| Sterility | Non-Sterile |
| Cuff | Beaded Cuff |
| Sizes | SM, MD, LG, XL, 2XL & 3XL |
| Accelerators | Yes |
| Chemo Tested | Yes |
| Glove Interior | Powder-Free, pH Balance |
| Glove Exterior | Full Textured |
| Length (inches) | 12 in |
| Palm Thickness (Avg) | 5.9 mil |
| Finger Thickness (Avg) | 8.0 mil |
| Tensile Strength (Mpa) | 35 / 30 Mpa |
| Ultimate Elongation (%) | 680 / 650 |
| Pin Hole (AQL) | 1.5 |



Quality Standards:

Manufactured in accordance with Quality System ISO 9001.
Exceeds and meets current ASTM Standards: - D6319, D5151, D1671 & D6978.
This product is made from Synthetic Nitrile Polymer and contains no natural latex rubber.
* When used in accordance to the guidelines set forth in the USP 800 chapter.
Tested chemotherapy drugs available upon request..

