




SCATTERGUARD

FLUOROSCOPIC PATIENT DRAPES

PROTECT YOURSELF FROM SCATTER RADIATION

During procedures requiring fluoroscopy, direct beam radiation contacts a patient and “scatters” towards the surgeon and other personnel. **“The scattered radiation from the patient comprises the main source of radiation dose to staff.” – International Atomic Energy Agency**

Protect yourself and your team with Scatterguard drapes. Tested in accordance with ASTM and IEC standards, ScatterGuard products provide excellent protection for fluoro-lab and operating suite personnel during interventional procedures. They are lightweight, lead free, conform to the patient’s body and are disposable through normal waste procedures.



PROTECTION PROPERTIES

LEAD EQUIVALENCE (LE)	ATTENUATION AT 90KVP
0.125mm	>= 75%
0.25mm	>= 90%
0.375mm	>= 95%



LEAD FREE



LATEX FREE




STERILE



DISPOSABLE



FLEXIBLE



CERTIFIED TO IEC AND ASTM STANDARDS



IEC 61331-1:2014



ASTM F3094-14

