

Vee-Shields

Introducing a NEW line of disposable patient eye shields.

Vanalay disposable eye shields are available in 3 protective options: **Laser, IPL and LED**. The comfort adhesive layer makes for easy application, adjustment and removal without damaging eyelashes. **Vee-Shields** are latex free and hypoallergenic.

› **Laser Vee-Shield** REF #: 816023 24 pair

Disposable **Laser Vee-Shields** provide patient protection from all medical lasers from 190nm-11,000nm. An added heat protective layer provides additional protection if a laser strike should occur. The comfort adhesive layer makes for easy application, adjustment and removal without damaging eyelashes. **Laser Vee-Shields** are latex free and hypoallergenic.



› **IPL Vee-Shield** REF #: 816009 50 pair

Disposable **IPL Vee-Shields** provide patient protection from the bright light derived from Intense Pulsed Light (IPL) devices. The comfort adhesive layer makes for easy application, adjustment and removal without damaging eyelashes. **IPL Vee-Shields** are latex free and hypoallergenic.



› **LED Vee-Shield** REF #: 816016 50 pair

Disposable **LED Vee-Shields** provide patient protection from the intense blue, yellow and red light of LED devices. The comfort adhesive layer makes for easy application, adjustment and removal without damaging eyelashes. **LED Vee-Shields** are latex free and hypoallergenic.



Manufactured BY:

Vanalay, LLC