

Creating A Better Care Experience For Women

KLEAN KATCH

Disposable Urine Specimen Cup Holder

Klean Katch is proud to be a new product that helps create a better healthcare experience for women. We are based and manufactured in the United States.



kleankatch.com

U.S. Patent No. 10,349,923

Creating A Better Care Experience For Women

KLEAN KATCH creates a better healthcare experience for women by **dignifying the collection** process of bodily materials (urine)

Problem

When a woman is asked to provide a urine sample, she is handed a specimen cup and asked to catch the urine midstream. This often results in a messy and unsanitary experience for the patient.

It is especially difficult for pregnant women, larger women, those with short arms or decreased flexibility to reach down and hold the specimen cup in the proper position.



Patients find this experience unsatisfactory. Nurses and staff may need additional time to clean up after the patient. The process of urine sample collection is often displeasing to all involved.

Solution **Products Created by Women**

KLEAN KATCH is a simple to use, low cost, disposable device that cradles the urine sample cup in it's base and has a convenient handle for women to hold on to. This leaves the hands away from the urine stream, making the experience more hygienic. Klean Katch transforms a distasteful experience into one that is satisfactory for the patient and staff.

- Creates positive feedback from women
- Eliminate mess for patients and staff
- Low cost solution for customer satisfaction
- Also great for pediatrics and elderly
- Optimal design for strength and function
- Minimal storage needed
- Easily disposable
- Limits contamination



Testimonials

I am pregnant and have a big belly. I have trouble reaching to give A sample. This really helped.
-Bev D, 26 yrs old.

This is the best idea I have seen in a long time. Keeping your hands clean is great. Nice to have a handle.
-Shawna C, 56 yrs old.

I hate giving urine samples, I seem to always make a mess. This was simple and clean, great job. Thanks.
-Alexis L, 30 yrs old.

This device made my urine collection so much easier. It has been a challenge in the past.
-Heather B, 49 yrs old.



KLEAN KATCH can customize the holder for ANY SPECIMEN CUP you use, please call us for details.



Survey Information



- Used the Product without experiencing urine on hands
- Say they were able to collect a full clean sample
- Say the experience is better using the Klean Katch system
- Would like doctors to carry the Klean Katch device
- Say they would recommend to a friend



Creating A Better Care Experience For Women

KLEAN KATCH

Disposable Urine Specimen Cup Holder

Single Holder Options

Item Number KK001-50ST

Holder Construction



One standard Klean Katch holder fits 4 different size specimen cups.



Creating A Better Care Experience For Women

KLEAN KATCH

Disposable Urine Specimen Cup Holder

Klean Katch Details

Item Number KK001-50ST



This is the Klean Katch Holder Fully Un-Folded and ready to be used by patient.



- **Folds flat for easy storage in drawer or on a shelf.**
- **Able to hold up to a 2.0625 inch cups.**
- **Designed and made in the USA.**
- **Printed using a 4 Color Process + Aqueous Coating.**
- **Made with 28lb cardboard.**
- **Shrink wrapped packages of 50.**