

The PhysioPRO



The **PhysioPRO Systems** are clinical grade shockwave therapy systems that induce an intracellular and extracellular reactions to assist in the healing and recovery of soft tissue.

Shockwave technology theory works by a deep penetrating focused wave that activate cells, driving the formation of new cells that replace damaged tissues. Shockwave therapy accelerates blood flow to the treated area and may induce vasodilation and angiogenesis restoring and improving musculoskeletal tissue.

- Clinical Grade 60-185 mJ (1/5 bar) / 1-22 Hz
- Deep Penetrating Focused Shockwave
- Pain Free Quick Treatments
- Preprogrammed Therapy Protocols
- Easy to Operate User Interface

The PhysioPRO: Portable Pulse Wave Therapy for Pain Management

Pre-Programmed Therapy Protocols for:

- Radial And Ulnar Epicondylitis
- Tendinitis Of The Shoulder
- Acute Calcific Tendinitis Of The Shoulder
- Post Muscular Injury
- Patellar Tendinopathy
- Achillodynia
- Plantar Fasciitis
- Heel Spurs
- Myofascial Trigger Point Therapy - Neck
- Myofascial Trigger Point Therapy - Back
- Trochanteric Bursitis
- Periostitis / Shin Splints
- Dupuytren's Disease
- Thumb Basal Joint Arthritis / Rhizarthritis

System Features & Benefits:

- 7.0" Color Touch Screen
- "Plug & Play" System With Touch Screen
- Multi-Preset Treatment Protocols
- User Friendly Interface
- Low Cost Of Ownership

The PhysioPRO TECHNICAL SPECIFICATIONS:

System Type	Main Unit with Therapy handle	Protocols	More than 25 illustrated preset treatment recommendations
Technology	Compressor free ballistic radial pulse therapy-system with electromagnetic generator	Controls	Color touch-screen for all software operations
Power Levels / Energy	60 to 185 mJ (equivalent to 1-5 bar)	Applicators	6/15/25mm
Frequency	1-22 Hz	Life time	Minimum warranty of 4,000,000 shocks per handpiece Minimum warranty of 150,000 shocks per applicator head
Modes	4 Burst modes (continuous / 4 / 8 / 12 pulses)	Dimensions	12 x 9.5 x 5 Inches (L/W/H)
Programs	7 preset programs, adjustable	Weight	2.07kg (controller unit only)
Package	Aluminum Carrying Case	Warranty	3 years for controller unit and 1 year for accessories

The PhysioPRO STANDARD ACCESSORIES:

ORDERING CODE	DESCRIPTION	QUANTITY / PC.	ORDERING CODE	DESCRIPTION	QUANTITY / PC.
OCE001010101Z	Therapy Handle	1	OCE00102040001	Power Line	1
OCE001020102	Foot Switch	1	OCE001010101Z-04	Handle Pedestal	1
OCE001010101Z-07	6mm Therapeutic Head	1	OCE00101010405	Stylus Pen	1
OCE001010101Z-08	15mm Therapeutic Head	1	OCE001010101Z-06	Silicone Cap	10
OCE001010101Z-09	25mm Therapeutic Head	1	OCE00103040101	Aluminum Alloy Suitcase	1

Florida Medical Sales

4320 A1A South, Suite 2, Saint Augustine, FL 32080

admin@floridamedicalsales.com | Tel: 888.269.1810