



TuffSkin Blu Latex Exam Gloves

12-inch Cuff and 18 mil Thickness

Powder-Free Non-Sterile Latex Examination Gloves

Ideal for chemotherapy compounding and other related chemotherapy applications. These 18-mil latex specialty glove is designed to offer exceptional added coverage and increased protection where there is potential exposure to sharps and chemicals.

Characteristics and Overview:

- Diamond textured gloves provide superior durability to handle your toughest applications, wet or dry.
- Delivers superior protection and are ideal for decontamination, chemotherapy compounding, cleaning and sterilizing instruments and meet the demands in law enforcement and industrial safety.
- Tested for puncture resistance and comfortable feel.
- Meets USP 800 compliant and D6978-05 Standards for use of chemotherapy gloves when handling hazardous drugs* .



Item No.	Description	Size	Packaging
MG12181	TuffSkin Blu Latex PF Exam Glove	Small	50 eaches/box 10 boxes/case
MG12182	TuffSkin Blu Latex PF Exam Glove	Medium	50 eaches/box 10 boxes/case
MG12183	TuffSkin Blu Latex PF Exam Glove	Large	50 eaches/box 10 boxes/case
MG12184	TuffSkin Blu Latex PF Exam Glove	X-Large	50 eaches/box 10 boxes/case
MG12185	TuffSkin Blu Latex PF Exam Glove	2X-Large	45 eaches/box 10 boxes/case

Glove Specifications:

Brand	TuffSkin Blu
Material	Latex
Color	Dark Blue
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	SM, MD, LG, XL and 2XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free, Double Chlorinated
Glove Exterior	Diamond Textured
Length (inches)	12 in
Palm Thickness (Avg)	15.0 mil
Finger Thickness (Avg)	18.2 mil
Tensile Strength (Mpa)	38 / 35 Mpa
Ultimate Elongation (%)	780 / 750
Pin Hole (AQL)	1.5

Quality Standards

Manufactured in accordance with Quality System ISO 9001.

Exceeds and meets current ASTM 3578 standard for critical defects.

Caution: No sensitizing threshold to natural rubber latex has been established. Any individual with hypersensitivity to natural rubber latex should consult a physician before using any latex product.

